



LeSONO

Wireless & Handheld Veterinary Ultrasound Scanner



Wireless design



Fast connection on Android, iOS and Windows



Full scanning mode: B, M, CF, PD, PW*1



Compatible with DICOM 3.0 System*2



2 hours battery life*3



IP67 waterproof tested

*1 & 2. Optional

*3. With setting FPS as "power saver"



LV128 VET Series



Small and Medium Animal Scanner

Frequency Bandwidth (MHz) 3.6 ~ 8.5

Max Depth (cm) 12.6

Field of View 100

Applications Canine, Feline, Porcine, Ovine, Bovine, Mouse, Rabbit

Also provide Linear, Convex and Phased array



Equine



Large Animal



Cardiac

Clinical Images



Canine




Feline



Porcine



Ovine



Equine

Available on Android & iOS and Windows



Learn more

LEL TEK

+886-2-2913-7577

info@leltek.com

6F.-3, No. 293, Sec. 1, Beixin Rd., Xindian Dist., New Taipei City 231, Taiwan (R.O.C.)



LeSONO

Wireless & Handheld Veterinary Ultrasound Scanner



	LV128M	LV128L	LV128C	LV128PA
Frequency Bandwidth	3.6 – 8.5 MHz	4.2 – 12.5 MHz	2.0 – 5.0 MHz	1.3 – 3.7 MHz
Max Depth	12.6 cm	12.6 cm	30.0 cm	30.0 cm
Field of View	100°	-	60°	90°
Applications	Can be customized in app			
	<ul style="list-style-type: none"> • Vet small • Vet medium • Small canine • Medium canine • Feline • Porcine • Ovine • Bovine • Mouse • Rabbit 	<ul style="list-style-type: none"> • Small canine • Medium canine • Large canine • Feline • Porcine • Ovine • Bovine • Equine • Mouse • Rabbit 	<ul style="list-style-type: none"> • Large canine • Porcine • Ovine • Bovine • Equine 	<ul style="list-style-type: none"> • Abdomen • Cardiac
Wi-Fi	802.11 a/b/g/n (2.4G, 5G)			
Modes	B-Mode, M-Mode, Color Doppler*, PW*, Power Doppler*			
DICOM	DICOM 3.0*			
Battery Working Time	Up to 2 hours			
Operation System	 	 	 	

(*Optional)

- This product, features, options, and transducers are not commercially available in all countries.
- Sales and Shipments are effective only after the approval by the regulatory affairs. Please contact your local sales representative for further details.
- 5GHz Wi-Fi connection is not available in all countries.



Learn more

LEL TEK

+886-2-2913-7577

info@leltek.com

6F.-3, No. 293, Sec. 1, Beixin Rd., Xindian Dist.,
New Taipei City 231, Taiwan (R.O.C.)