

# Dual-Probe Veterinary Ultrasound Wireless Handheld Scanner



Wireless design



Fast connection on Android, iOS, Windows



Full scanning mode: B, M, CF, PD, PW



Compatible with DICOM 3.0



2.5 hours battery life



IP68 WaterProof



Canine



Feline



Porcine



Ovine

Access 10+ presets for quick and easy start or create your own custom presets.



Small Animal



Medium Animal



Small Canine



Medium Canine



Large Canine



Feline



Porcine



Ovine



Bovine



Equine



Hi-Res



	LXV192LC		LXV128LP		LXV128LC	
Transducer	L 12-5	C5-2	L 12-5	P4-2	L 12-5	C5-2
Frequency	4.2 - 12.5 MHz	2.0 - 5.0 MHz	4.2 - 12.5 MHz	1.3 - 3.7 MHz	4.2 - 12.5 MHz	2.0 - 5.0 MHz
Bandwidth Max	12.6 cm	30 cm	12.6 cm	30 cm	12.6 cm	30 cm
Depth Field of View	-	60°	-	90°	-	60°
Number of Element	192	192	128	64	128	128

Indications for Use and Applications (Customizable in App)

Small Canine	✓	-	✓	-	✓	-
Medium Canine	✓	-	✓	-	✓	-
Large Canine	✓	✓	✓	-	✓	✓
Feline	✓	-	✓	-	✓	-
Porcine	✓	✓	✓	-	✓	✓
Porcine Lard	-	✓	-	-	-	✓
Ovine	✓	✓	✓	-	✓	✓
Bovine	✓	✓	✓	-	✓	✓
Equine	✓	✓	✓	-	✓	✓
Mouse	✓	-	✓	-	✓	-
Rabbit	✓	-	✓	-	✓	-
VetMiddle	-	✓	-	✓	-	✓
Cardiac	-	✓	-	✓	-	✓

Wi-Fi	802.11 a/b/g/n (2.4GHz, 5GHz)					
Modes	B-Mode, M-Mode, Color Doppler*, PW*, Power Doppler*					
DICOM	DICOM 3.0*					
Wireless Charging	Support 5V and 9V, battery with 3000mAh					
Usage Hours	Up to 2.5 hours					



(\*optional)

- This product, features, options, and transducers are not commercially available in all countries.
- Sales and Shipments are effective only after the approval by the regulatory affairs.  
Please contact your local sales representative for further details.
- 5GHz Wifi connection is not available in all countries.



Learn more

**LEL TEK**

+886-2-2913-7577

info@leltek.com

6F.-3, No. 293, Sec. 1, Beixin Rd., Xindian Dist.,  
New Taipei City 231, Taiwan (R.O.C.)